

HW 1

Muhammad Elarbi

Summer 1 2024

Class Information

Class: EECE 4649: Biomedical Imaging

Instructor: Charles A. DiMarzio

Semester: Summer 1 2024

Name: Muhammad Elarbi

Solution(s)

Problem 1

(a)

Beer's Law:

$$P(z) = P(0) \cdot e^{-\mu_a \cdot z}$$

where:

- $P(z)$: The power of the light (or X-ray) after traveling a distance z through the medium. This is the power detected after attenuation.
- $P(0)$: Initial power (before entering the medium).
- μ_a : Absorption coefficient
- z : Distance to be solved. This is the distance that the light travels through the medium. In our context, this is the depth at which the mirror is placed in the fluid.

What We Are Solving For:

We are solving for z , which represents the maximum depth (distance we are solving for) at which the mirror can be placed below the surface of the fluid such that the reflected light can still be detected.

Given Data From The Problem:

- Initial power, $P(0)$: $0.1 \text{ mW} = 0.1 \times 10^{-3} \text{ W}$
- Minimum detectable power (from Detector), $P(2z)$: $1 \text{ nW} = 1 \times 10^{-9} \text{ W}$
- Wavelengths: 650 nm and 450 nm

Beer's Law for Round Trip:

Since the problem states the light travels to the mirror and back, **the distance is $2z$** . This changes the Beer's Law equation we are working with to be:

$$P(2z) = P(0) \cdot e^{-\mu_a \cdot 2z}$$

where:

- $P(2z)$: Detected power after traveling to the mirror and back.

Rearranging Beer's Law to Solve for z :

By rearranging Beer's Law and using the natural logarithm, we can solve for the maximum depth z at which the mirror can be placed. All we need to do is perform some algebraic manipulation and isolate z on one side of the equation:

1. Start with the equation:

$$P(2z) = P(0) \cdot e^{-\mu_a \cdot 2z}$$

2. We can substitute in the given values for $P(2z)$ and $P(0)$ to make things easier:

$$1 \times 10^{-9} = 0.1 \times 10^{-3} \cdot e^{-\mu_a \cdot 2z}$$

3. Divide both sides by $P(0)$:

$$\frac{1 \times 10^{-9}}{0.1 \times 10^{-3}} = e^{-\mu_a \cdot 2z}$$

$$10^{-6} = e^{-\mu_a \cdot 2z}$$

4. Take the natural logarithm (\ln) of both sides to solve for the exponent:

$$\ln(10^{-6}) = \ln(e^{-\mu_a \cdot 2z})$$

$$-6 \ln(10) = -\mu_a \cdot 2z$$

5. Solve for z :

$$-6 \cdot 2.3026 = -\mu_a \cdot 2z$$

$$6 \cdot 2.3026 = \mu_a \cdot 2z$$

$$13.8156 = \mu_a \cdot 2z$$

$$2z = \frac{13.8156}{\mu_a}$$

$$z = \frac{13.8156}{2 \cdot \mu_a}$$

Now, we can use the above equation to solve for the z values in both water and blood at both wavelengths.

Calculation for Water (H₂O) at 650 nm:**Absorption coefficient, μ_a (Estimate visually from the graph):**

- For water at 650 nm: $\mu_a \approx 10^{-2.5} \text{ cm}^{-1} \approx 3.162 \times 10^{-3} \text{ cm}^{-1}$

Calculate z :

$$z = \frac{13.8156}{2 \cdot \mu_a}$$

$$z = \frac{13.8156}{2 \cdot 3.162 \times 10^{-3}}$$

$$z = \frac{13.8156}{0.006324}$$

$$z \approx 2184.8 \text{ cm} \approx 21.85 \text{ m}$$

Calculation for Water (H₂O) at 450 nm:**Absorption coefficient, μ_a (Estimate visually from the graph):**

- For water at 450 nm: $\mu_a \approx 10^{-1.8} \text{ cm}^{-1} \approx 1.584 \times 10^{-2} \text{ cm}^{-1}$

Calculate z :

$$z = \frac{13.8156}{2 \cdot \mu_a}$$

$$z = \frac{13.8156}{2 \cdot 1.584 \times 10^{-2}}$$

$$z = \frac{13.8156}{0.031697}$$

$$z \approx 364 \text{ cm} = 3.64 \text{ m}$$

Note - How to calculate Absorption Coefficients for Blood:

The assignment says to assume equal concentrations of Oxy-Hemoglobin (HbO₂) and Deoxy-Hemoglobin (Hb). Therefore, we can calculate the average absorption coefficient for both wavelengths by **averaging the values for Oxy- and Deoxy-Hemoglobin from the figure.**

Calculation for Blood at 650 nm:**Absorption coefficient, μ_a (Averaged):**

$$\begin{aligned}\mu_a (\text{HbO}_2) &\approx 10^3 \text{ cm}^{-1} \\ \mu_a (\text{Hb}) &\approx 10^2.5 \text{ cm}^{-1} \approx 3.16 \times 10^2 \text{ cm}^{-1} \\ \mu_a (\text{average}) &\approx \frac{10^3 + 3.16 \times 10^2}{2} = 5.58 \times 10^2 \text{ cm}^{-1}\end{aligned}$$

Calculate z :

$$\begin{aligned}z &= \frac{13.8156}{2 \cdot \mu_a(\text{average})} \\ z &= \frac{13.8156}{2 \cdot 5.58 \times 10^2} \\ z &= \frac{13.8156}{1116.0} \\ z &\approx 0.0123 \text{ cm} = 0.123 \text{ mm}\end{aligned}$$

Calculation for Blood at 450 nm:**Absorption coefficient, μ_a (Averaged):**

$$\begin{aligned}\mu_a (\text{HbO}_2) &\approx 10^3.5 \text{ cm}^{-1} \approx 3.16 \times 10^3 \text{ cm}^{-1} \\ \mu_a (\text{Hb}) &\approx 10^2.8 \text{ cm}^{-1} \approx 6.31 \times 10^2 \text{ cm}^{-1} \\ \mu_a (\text{average}) &\approx \frac{3.16 \times 10^3 + 6.31 \times 10^2}{2} = 1.89 \times 10^3 \text{ cm}^{-1}\end{aligned}$$

Calculate z :

$$\begin{aligned}z &= \frac{13.8156}{2 \cdot \mu_a(\text{average})} \\ z &= \frac{13.8156}{2 \cdot 1.89 \times 10^3} \\ z &= \frac{13.8156}{3780.0} \\ z &\approx 0.0036 \text{ cm} = 0.036 \text{ mm}\end{aligned}$$

Summary:

- **Water at 650 nm:** $z \approx 21.85$ m
- **Water at 450 nm:** $z \approx 3.64$ m
- **Blood at 650 nm:** $z \approx 0.123$ mm
- **Blood at 450 nm:** $z \approx 0.036$ mm

(b)

These results show that red light (650 nm) penetrates deeper than blue light (450 nm) in both water and blood. This is because red light is less absorbed by these media compared to blue light.

Problem 2

(a)

Given Data:

- Wavelength of X-ray beams, $\lambda = 0.1 \text{ nm}$
- Initial intensity of X-ray beams, I_0
- Phantom thickness, 10 cm
- Diameter of bone, 3 cm
- Mass Absorption Coefficient for soft tissue, $\mu/\rho = 5.38 \text{ cm}^2/\text{g}$
- Mass Absorption Coefficient for bone, $\mu/\rho = 28.51 \text{ cm}^2/\text{g}$
- Density of soft tissue, $\rho_{\text{tissue}} = 1.06 \text{ g/cm}^3$
- Density of bone, $\rho_{\text{bone}} = 1.92 \text{ g/cm}^3$

Calculation of Attenuation for I_1 and I_2 :

For beam I_1 :

The beam I_1 passes through only the soft tissue.

1. Path Length through Soft Tissue:

$$d_{\text{tissue}} = 10 \text{ cm}$$

2. Attenuation Coefficient for Soft Tissue:

$$\mu_{\text{tissue}} = \left(\frac{\mu}{\rho} \right)_{\text{tissue}} \times \rho_{\text{tissue}} = 5.38 \text{ cm}^2/\text{g} \times 1.06 \text{ g/cm}^3 = 5.7028 \text{ cm}^{-1}$$

3. Total Attenuation for I_1 :

$$I_1 = I_0 e^{-\mu_{\text{tissue}} d_{\text{tissue}}} = I_0 e^{-5.7028 \times 10}$$

$$I_1 = I_0 e^{-57.028}$$

For beam I_2 :

The beam I_2 passes through both the bone and the soft tissue.

1. Path Length through Bone:

$$d_{\text{bone}} = 3 \text{ cm}$$

2. Path Length through Soft Tissue (excluding bone):

$$d_{\text{tissue}} = 10 \text{ cm} - 3 \text{ cm} = 7 \text{ cm}$$

3. Attenuation Coefficient for Bone:

$$\mu_{\text{bone}} = \left(\frac{\mu}{\rho} \right)_{\text{bone}} \times \rho_{\text{bone}} = 28.51 \text{ cm}^2/\text{g} \times 1.92 \text{ g/cm}^3 = 54.7392 \text{ cm}^{-1}$$

4. Total Attenuation for I_2 :

$$I_2 = I_0 e^{-(\mu_{\text{bone}} d_{\text{bone}} + \mu_{\text{tissue}} d_{\text{tissue}})}$$

$$I_2 = I_0 e^{-(54.7392 \times 3 + 5.7028 \times 7)}$$

$$I_2 = I_0 e^{-(164.2176 + 39.9196)}$$

$$I_2 = I_0 e^{-204.1372}$$

Summary:

• **Total Attenuation for I_1 :**

$$I_1 = I_0 e^{-57.028}$$

• **Total Attenuation for I_2 :**

$$I_2 = I_0 e^{-204.1372}$$

(b)

Synopsis/Logic: The process begins by defining key parameters such as the thickness of the soft tissue, the diameter of the bone, and the attenuation coefficients for both materials. The code then creates a series of positions along the length of the tissue sample. For each position, the code checks if the position falls within the region containing the bone. If it does, the code calculates the attenuation by considering both the bone and the surrounding soft tissue. If the position falls outside the bone region, only the soft tissue's attenuation is considered. The attenuation is calculated using Beer's Law, represented by the exponential decay formula $I = I_0 e^{-\mu x}$, where μ is the attenuation coefficient and x is the thickness of the

material. Finally, the results are visualized using a grayscale plot, where the intensity of the color indicates the level of attenuation, with darker regions representing higher attenuation due to the presence of bone.

Code:

```
% X-ray Attenuation through Tissue with Bone
% This script calculates and visualizes the attenuation of X-
% rays as they pass through
% a 10 cm thick tissue sample containing a 3 cm diameter bone.
%
% Author: Muhammad Elarbi
% Email: elarbi.m@northeastern.edu
% Date Created: May 13, 2024
% Last Modified: May 14, 2024
% Course: EECE 4649: Biomedical Imaging
% Instructor: Charles A. DiMarzio
%
% Usage:
% This script models the attenuation of X-ray beams through a
% tissue sample with an embedded bone.
% Run the script in MATLAB to view the output in the Command
% Window.

% Define parameters
tissue_thickness = 10; % cm
bone_diameter = 3; % cm
soft_tissue_mu = 5.7028; % cm-1
bone_mu = 54.7392; % cm-1
I0 = 1; % Initial intensity (assumed to be 1 for simplicity)

% Create x positions from 0 to 10 cm
x = linspace(0, tissue_thickness, 1000);


% Initialize attenuation array
attenuation = zeros(size(x));

% Calculate attenuation for each position
for i = 1:length(x)
    if x(i) >= (tissue_thickness - bone_diameter) / 2 && x(i)
        <= (tissue_thickness + bone_diameter) / 2
        % Passes through bone
        attenuation(i) = I0 * exp(-bone_mu * bone_diameter -
            soft_tissue_mu * (tissue_thickness - bone_diameter))
        ;
    end
end
```

```
    else
        % Passes through soft tissue only
        attenuation(i) = I0 * exp(-soft_tissue_mu *
            tissue_thickness);
    end
end

% Plot the result
figure;
imagesc(x, [0 1], attenuation);
colormap('gray');
colorbar;
xlabel('Position (cm)');
ylabel('Intensity');
title('X-ray Attenuation through Tissue with Bone');
```

The figure below shows the result of the MATLAB code, visualizing the attenuation of X-rays through the tissue with an embedded bone.



m_2024_summer1_eece4649_elarbi_hw1_q2_b.pdf

Figure 1: X-ray Attenuation through Tissue with Bone

Problem 3

(a)

(b)

Problem 4

(a)

(b)

(c)

(d)